

SERVICE

OPPORTUNITIES

Flannel Diaper Sets

A matching set of flannel diapers is included in each Caring for Coy care package. They are soft against delicate, sometimes underdeveloped skin, and are used to offer dignity and modesty for baby for photographic memories. Angel babies range in size from micro-preemie to full-term weight, so a variety of sizes is needed.

Please use baby flannel for all diapers, as the material is more flexible and does not harm baby's fragile skin.




General Instructions:

1. Cut baby flannel in sets of two, using attached patterns.
2. With right sides together, stitch with 1/4" guide around entire piece, leaving 1" unstitched for turning.
3. Trim curves and turn right-side-out.
4. Press.
5. Top-stitch all around piece to finish unstitched area.
6. Fold center upward to meet back edge, both sides in toward center and secure with appropriate-sized safety pin.
7. Create identically-matched sets, so one can be placed lovingly on baby and one can be kept by the family in a memory box.
8. Tie sets together with soft ribbon or bag in pairs in mini Ziplok bags.




Volunteer Services
UINTAH BASIN HEALTHCARE


Flannel Diaper Patterns Small & Medium



Volunteer Services
UINTAH BASIN HEALTHCARE



Flannel Diaper Patterns
Large



Volunteer Services
UINTAH BASIN HEALTHCARE



Brody's Fund
UINTAH BASIN HEALTHCARE

